

# VRA Centring Toy

September 2021

## Glo-worm

- Glo worm face lights up when light is activated from the VRA console



(In stock design and colour choices vary)

## Optional Acrylic case for Centring toy

- Conceals toy until light is activated, made of dark translucent acrylic
- Acrylic case is a sanitizable surface



Dimensions ( approximate): H 28cm, W 20cm, D 11cm

### Care of Surface:

Use a wet microfibre cloth to remove dust and finger prints

Cleaning wipes with at least 70% alcohol, hydrogen peroxide, bleach, isopropyl alcohol or ethyl alcohol, or other acrylic approved disinfectants.

Never use window cleaners or other disinfectants that contain ammonia, such as Windex

(Clean with high pressure as Acrylic scratches)



electro-medical instrument co.

1-2359 Royal Windsor Drive, Mississauga, ON L5J 4S9

800-263-6430, 905-822-3188

[emi-canada.com](http://emi-canada.com)

Your source for audiology equipment, supplies and services, since 1972